

B.safe waste cap electronic fill-level control



To prevent your waste canister from overflowing, connect these [b.safe](#) waste caps with electronic fill-level control to the alarm box. When the 'canister full' level is reached, you will receive an optical and acoustic warning. The alarm is triggered before your canister is completely full. This leaves a small safety buffer in the canister.

- Prevent your waste canisters from overflowing
- Optical and acoustic warning
- 2 alarmbox models available
- Accessories like fittings and hose connectors guarantee you waste lines are connected leak tight
- Special couplings for corrugated flexible tubes
- Exhaust filters for different sizes of waste containers adsorb hazardous vapours

Fittings, hose connectors and corrugated hose couplings allow you to configure the system according to your requirements. Different thread sizes guarantee that you will find the right cap for your canister.

The way of function is quite simple: connect your [b.safe](#) waste cap with the [b.safe](#) display by means of a suitable cable. With increasing level, the cap integrated floater rises. As soon as the canister is full, a contact is closed between the cap and the floater and an appropriate

signal is transferred to the [b.safe](#) display. You can easily recognize which container has to be exchanged through visual and acoustic signs on the [b.safe](#) display.

The alert already rings before the container is completely filled to have a little safety buffer. Even waste from the sample preparation can be directly poured into the waste container via a permanently mounted funnel.

A special activated carbon in the [b.safe](#) exhaust filters absorbs all solvent vapours in order to have only harmless air coming into the environment. Practical service life indicators guarantee full control on the permitted service life.



Plot No 39 1st Floor, Lavanya Aracade, Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.
Mob : +91 99082 46685 Tel : 040 - 35172288/89