



RS-SC-2 Automatic Seed Counting Machine

INTRODUCTION:

Automatic seed Counter can automatic count many kinds of grains by setting number. In agriculture it is mainly used in grain count, exp. Rice, wheat, sorghum, corn, and vegetables seed and so on.

FEATURES

1. Two modes: Free count and bulk count
2. The counting speed is adjustable
3. High accuracy: Sensitivity levels for different sizes of particles
4. Automatic deceleration: when the counting number is close to the setting number, the counting speed slows down automatically
5. Automatic motor Protection: the motor will stop working once it detects there are no particles on the plate for 60seconds
6. User-friendly interface: LCD screen, mode selection, counting time, level, and number of particles displayed clearly
7. Wide application: suitable for particles includes but is not limited to crop seeds, vegetable seeds, feed, diamond, and tablets



TECHNICAL PARAMETER:

Count Range	0.7 - 12.0 mm
Count Accuracy	± 2%
Count Speed	≥1000/3 min (Sesame)
Counting Capacity	1~99999
Bulk Counting	1~99999
Noise	<70 db
Power	AC110V/60 Hz or AC220V/50Hz
Continuous Working Time	≥ 5 Hr
Size of Machine	285*210*176 mm
Operation Environment	-10° C-50° C, <85% RH



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.