

DIGITAL COLONY COUNTER



Specifications

Digital colony counter is designed for quick and accurate counting of bacterial and mould colonies in petri dishes. Our digital colony counter is a feature packed ft easy to use, tins is an indispensable bench top tool for the busy microbiologist. It is designed for rapid and accurate counting of bacterial and mould colonies.

- 100 mm Dia Magnifying Glass
- Audible Confirmation of Each Count
- Uniform Glair Free illumination
- Digital Readout from 0-9999
- Auto Marker Pen
- Wolffhuegel Glass Grid with Focusing Facility

Display	4 Digit, 9999 Max. Count
Dish Size	110 mm
Magnification	x 1.7
Dimension	L 274 X B 320 X H 167 mm (Approx)
Weight	3 Kgs. (Approx)
Packing	23cm x 28cm x20cm Corrugated Carton
Power	230V – 10% AC, 50 Hz, 40W
Accessories	Marking Pen – 1 No.
Accessories	Magnifier Lens – 1 No.

