

TABLET DISSOLUTION TEST APPARATUS



Technical Specifications:

- 8 Nos. Easy to attach and detach `Press Fit` wobble free Stirrer shafts.
- 8 Nos. simple attachable or detachable Paddles Assembly (USP II)
- 8 Nos. Baskets with simple attachable or detachable Assembly(USPI)
- 8 Nos. Clear Dissolution Glass Vessels.
- Auto Lifting Arrangement
- Water Bath with non-corrosive Tank plate.
- Large Graphic LCD (240 X 64) to indicate set parameters like: RPM, Vessel
- Temperature, time interval, test time, etc.
- Wobble and vibration free instrument.
- Precise individual vessel cantering system.
- Facility to review set parameters.
- Non-volatile memory.
- Memory to store 48 test records.
- RPM: 25 TO 160.
- Temperature Range: 30 .0°C – 45. 0°C with accuracy $\pm 0.1^{\circ}\text{C}$.
- Tablet dropping system of SS provided to insert tablets simultaneously.
- Lids provided for minimal evaporation
- IQ, OQ documents will be provided with the instrument without any charge.

