



Dissolved Oxygen Meter

RS-PSI 1601

DESCRIPTION:

Digital Dissolved Oxygen Meter RS- 801 & 811 are available in benchtop models. These instruments are extremely useful for agriculture, soil analysis, fisheries, water quality control in boiler feed water, water works, breweries, water purification plants, chemical and pharmaceutical industries etc.

FEATURES:

1. Highly Stable & Accurate
2. 3½ Digit LCD Display
3. DO & Temperature Measurements



TECHNICAL SPECIFICATIONS:

Model	RS 801	RS 811
Range	DO: 0-20 ppm Temp: 0-60° C	DO: 0-20 ppm -
Resolution	DO: 0.1 ppm Temp: 0.1° C	DO: 0.1 ppm
Accuracy	DO: 0.1 ppm Temp: 0.2° C	DO: 0.1 ppm -
Stirrer	Available	-
Temperature Compensation	0 - 50 ° C Manual	
Calibration	Manual	
Display	3½ Digit Seven Segment LED	
Power	230V AC ± 10%, 50 Hz	
Dimensions	275X135X75 mm (IXbXh) approx.	
Weight	3.5 Kgs approx.	
Accessories	Gold/Silver Amperometric DO probe, Set of DO membranes, Temperature probe (Model-801), Stirrer and Magnetic Capsule (Model-801), Operation Manual and Dust Cover	

Note: Specifications are subject to change due to continuous improvement.



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.