

HPLC COLUMN WASHING PUMP SA100



Specifications:

- Flow Rate : 0.1 to 10.00 mL/min for 10 mL/min head
- Pressure : 0 to 6,000 psi
- Flow Accuracy : +- 5% for a Flow rate of 0.1 to 10.000 mL/min, with 100% Water @1000 psi
- Dimensions : 5.5" High x 10.375" Wide x 17.5" deep
- Weight : 24 lb
- Power : 100-120 VAC, 50-60 Hz; or 220-240 VAC, 50-60 Hz
- Features : 4 Line LCD Display, Splash Proof Keypad, Prime Purge Valve, Auto Stop Timer, Service Mode
- Optional : Pulse Damper
- Nos of Column : Upto 4 Column Can be Washed Simultaneously
- Remote Inputs : RS-232

