



RSolv[®]

LIFE SCIENCES

LAB CLEAN

Best Solution For Your Laboratory Glassware



Lab Clean generally used in Research centers, Path labs, Diagnostic centers, Hospital and many more.

Ingredients :

LABSA/DBSA -2%, SODA ASH -2%, STP -20%,
SLES -5% Aqua Solution -81%

Method of Use :

Pour 70 ml of Lab Clean in a bucket and drop all glassware/plasticware. Take out after 15min and wash with regular water



| | |
|-------------|---------------------------------|
| Mfg. Date | : Aug 2023 |
| Expiry Date | : One year from the date of mfg |
| Batch No. | : R-008 |
| ML No. | : RAK-2721 |
| M.R.P. | : 1500.00 |

Distributed By:

RIKI GLOBAL TRADING PVT LTD

Ph.No : 040 - 42217787, 9908246695

Email Id : info@rikiglobal.com

Website : www.rikiglobal.com