

Bagmixer

Bagmixer 400 w



Product Specification:

- Brand Interscience
- Model Number RS-BagMixer 400 W
- Capacity 50-400 mL
- Material Stainless Steel
- Blending speed Fixed speed 8 strokes/s
- Dimensions 26 x 39 x 29 cm
- Weight 16.5 kg
- Power 110-220 V 50-60Hz
- Blending time 30,60,90,120,150,180,210 s
- Compact and ergonomic Yes
- Leak free Q tight system Yes
- Window door Yes
- Security drip tray Yes
- Door 270 deg opening door
- Door weight 950 g

Product Description:

Avoiding all risks of cross-contamination, RS-BagMixer 400W is an easy-to-use a powerful lab blender. Adapted to all kind of applications and with a gurantee of optimal bacterial extraction, it is a great tool for lab analyses.

1. QUICK HOMOGENIZING OF SAMPLES

- Optimal bacterial extraction
- Only 60 seconds to homogenize a sample

2. EASY CLEANING

- Fully opening door (270°)
- Security drip tray
- All stainless steel

3. USER-FRIENDLY

- Compact
- Adjustable blending time (30 -210 secs. or open-running)
- Q-Tight® closing: no more leaks upwards

Plot No 39 1st Floor Lavanya Arcade Jayabheri Enclave Gachibowli Hyderabad-500032.Telangana
Mob : +91 99082 46685 Tel : 040 - 35172288/89