



Tablet Dissolution Test Apparatus

TECHNICAL SPECIFICATIONS:

1. 8 Nos. Easy to attach and detach 'Press Fit' wobble free Stirrer shafts
2. 8 Nos. simple attachable or detachable Paddles Assembly (USP II)
3. 8 Nos. Baskets with simple attachable or detachable Assembly(USPI)
4. 8 Nos. Clear Dissolution Glass Vessels.
5. Auto Lifting Arrangement
6. Water Bath with non-corrosive Tank plate.
7. Large Graphic LCD (240 X 64) to indicate set parameters like: RPM, Vessel
8. Temperature, time interval, test time, etc.
9. Wobble and vibration free instrument.
10. Precise individual vessel cantering system.
11. Facility to review set parameters.
12. Non-volatile memory.
13. Memory to store 48 test records.
14. RPM: 25 TO 160
15. Temperature Range: 30 .0°C – 45. 0°C with accuracy $\pm 0.1^{\circ}\text{C}$
16. Tablet dropping system of SS provided to insert tablets simultaneously.
17. Lids provided for minimal evaporation
18. IQ, OQ documents will be provided with the instrument without any charge.



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.