

RSolv Viscometer is a scientific instrument used to measure the viscosity of a fluid. Viscosity is a measure of a fluid's resistance to deformation or flow, and it describes how easily a fluid can flow. In simple terms, it is a measure of a fluid's thickness or stickiness.



RSolv Viscometer

Features

- Full color, touch screen display
- Accurate and reliable operation
- Direct readout of all measurement parameters
(Viscosity, Temperature, % Torque, Spindle/Speed)
- Auto Range Showing
- Time Function for measurement
- Sound alarm at under 20% Torque
- Linear calibration

