

# VACCUM FILTRATION PEDESTAL

HANDS-FREE MEDIA FILTERING MADE EASY AND AFFORDABLE



Use in Conjunction with the available separately Autofil® filtration system and/or the Cradle Ring. These two filtering options offer a permanent vacuum hose connection and minimize set-up time. The vacuum hose connects to the back of the pedestal or cradle ring and enables adjustable flow control using the on/off switch. This patented docking system is made from aluminum.

## PATENTED DOCKING HARDWARE:

- Allows for hands-free filtration
- Provides on/off vacuum flow control
- Permanent hose connection to minimize set-up time
- Works in conjunction with Autofil® system in 15ml, 50ml, 250ml, 500ml, and 1000ml

## SPECIFICATIONS:

- Color: Black
- Made in: USA
- Material: Aluminum
- Size/Dimensions;
- Overall height: 8"
- Base diameter: 7"
- Base thickness: 1.5"
- Clearance - base to dock: 5.75"

