

Digital Overhead Stirrer

Product Features:

- Independent LED display for RPM & Timer 99:59 (HH:MM)
- User Friendly
- Brushless DC motor
- Liquid cannot reach internal components due to the front safeguard
- Audible alarm for error conditions
- Effective stirring provides clockwise & anti-clockwise stirring modes
- Prominently used in chemical synthesis, pharmaceuticals, physical & chemical analysis, petrochemical, cosmetics, health products, foods, biotechnology & so on



Technical Specification:

Model No.	6100.RSolv	6050.RSolv
T-On Time	: 5 – 999 Sec	: 5 – 999 Sec
T-Off Time	: 0, 5 – 999 Sec	: 0, 5 – 999 Sec
Speed Range	: 50 – 3000RPM	: 50-1500RPM
Wifi Control	: NO	: NO
Temperature Measurement	: NO	: NO
Max. Viscosity	: 6000 mPass	: 10000 mPass
Max. Torque	: 65N.CM	: 80N.CM
Time Range	: 1min – 99:59min	: 1min – 99:59min
Max. Stir Capacity (H2O)	: 10-40L	: 40-60L
Stirring Shaft Extension Feature	: Yes	: Yes
Stirring Modes	: Yes	: Yes
Motor Type	: BLDC	: BLDC
Operating Voltage	: 24 VDC	: 24 VDC
Supply Voltage	: 100 – 240Vac, 50/60Hz	: 100 – 240Vac, 50/60Hz
Power	: 130	: 130
Dimension	: 90x240x220mm	: 90x240x220mm
Net Weight	: 12 Kg	: 12 kg

