

## Solvent Waste System



### SPECIFICATIONS:

1. Carbon Exhaust Filler traps' solvent vapors
2. Tube Fittings secure waste tubes to container in order re- create a closed-loop system. preventing spills and vapors from leaking
3. 90 degree Adapter included for maximum storage
4. Versacap technology eliminates tube tangling when interchanging open and closed cap
5. EPDM gasket for high chemical compatibility against aggressive organic solvents
6. Includes plugs for unused for ports

Plot No 39, 1st Floor, Lavanya Arcade,  
Jayabheri Enclave, Gachibowli,  
Hyderabad-500032.Telangana.