

DIGITAL GLASS BEAD STERILIZER



Features:

- Portable size
- Rapid warm-up facility
- High efficiency and highly accurate temperature
- Digital micro-controlled based
- Ready to use signal with buzzer
- Auto sleep mode

Specifications:

- Dimension (mm): 200 x 15 x 36
- well Size mm 50x58
- Wattage:250w
- Temperature Controller Digital Temp indicator with Micro Controller based controller
- Weight incl Glass Beads: 21 Kg
- Input Power: 220-240 V, 50 Hz

Plot No 39 1st Floor Lavanya Arcade Jayabheri Enclave Gachibowli Hyderabad-500032.Telangana
Mob : +91 99082 46685 Tel : 040 - 35172288/89