

Digital Autoclave

SPECIFICATIONS:

1. Single wall unit constructed out of stainless steel
2. **Safety devices:** Safety valve, weight valve
3. Micro Processor auto-tune PID controller with timer
4. Standard features: steam release valve, pressure gauge, stand for bin, silicone/neoprene gasket
5. **Heating element:** 2000 watts

Size in cms
30X30
30X50
30X60
30X70



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.