

# Vertical Autoclave



	<b>DIAX HEIGHT</b>	<b>Ltr.</b>	<b>LOAD</b>
1.	30 cm x 50 cm	35	2KW
2.	35 cm x 55 cm	52	3 KW
3.	40 cm x 60 cm	80	4 KW
4.	45 cm x 60 cm	98	4KW
5.	45 cm x 75 cm	120	5KW
6.	55 cm x 75 cm	152	6 KW
7.	55 cm x 90 cm	200	6 KW

## TEMPERATURE RANGE:

- 121 C to 125 C
- Temperature Resolution 0.1 C
- Temperature Accuracy  $\pm 0.5$  C
- Temperature Sensor
- PT100 RTD type Class 'A'
- Temperature Control
- Low Water Level Alarm
- Microprocessor based PLC digital display.
- Auto purging and auto release
- Temperature Display
- LED / LCD back lit display.

## PRESSURE RANGE:

- 15 to 22 PSI, factory set at 15PSI
- Pressure Display
- Dial Gauge
- Pressure resolution
- Electrical
- 230 V / 440 V / 50Hz

Plot No 39 1st Floor Lavanya Arcade Jayabheri Enclave Gachibowli Hyderabad-500032.Telangana  
Mob : +91 99082 46685 Tel : 040 - 35172288/89