



# RSolv Washing Pump

## SPECIFICATIONS:

Flow Rate	0.1 to 10.00 mL/min for 10 mL/min head
Pressure	0 to 6,000 psi
Flow Accuracy	± 5% for a Flow rate of 0.1 to 10.000 mL/min, with 100% Water @1000 psi
Dimensions	5.5" High x 10.375" Wide x 17.5" deep
Weight	24 lb
Power	100-120 VAC, 50-60 Hz; or 220-240 VAC, 50-60 Hz
Features	4 Line LCD Display, Splash Proof Keypad Prime Purge Valve, Auto Stop Timer, Service Mode
Optional	Pulse Damper
Nos of Column	Upto 4 Column Can be Washed Simultaneously
Remote Inputs	RS-232



Plot No 39, 1st Floor, Lavanya Arcade,  
Jayabheri Enclave, Gachibowli,  
Hyderabad-500032.Telangana.