



Vacuum Filtration Pedestal

Hands-Free Media Filtering
Made Easy and Affordable

FEATURES:

Use in Conjunction with the available separately Autofil® filtration system and/or the Cradle Ring. These two filtering options offer a permanent vacuum hose connection and minimize set-up time. The vacuum hose connects to the back of the pedestal or cradle ring and enables adjustable flow control using the on/off switch. This patented docking system is made from aluminum.

PATENTED DOCKING HARDWARE:

1. Allows for hands-free filtration
2. Provides on/off vacuum flow control
3. Permanent hose connection to minimize set-up time
4. Works in conjunction with Autofil® system in 15ml, 50ml, 250ml, 500ml, and 1000ml



SEPCIFICATIONS:

Color	Black
Made in	USA
Material	Aluminum
Size/Dimensions	Overall height: 8" Base diameter: 7" Base thickness: 1.5"
Clearance - base to dock	5.75"



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.