



RIKI GLOBAL
performance through innovation

KJELDAHL DIGESTER



Features:

- Soft touch tactile membrane keyboard with LED indication

Specifications:

- Heating Block Solid: Aluminum
- No. of samples: 8 nos
- Glass Flask: 250 ml
- Temperature : 10° C above ambient to 400° C over temp
- Sample Size: 0.25 - 5 gms

Protection Available Temp: Resolution 0.1° C sensor PT 100

Timer: User definable up to 99 Hrs 59 Min.

Display: 7-segment LED keyboard soft touch membrane

Power: 230V Ac, 50 Hz wattage 1500 Watts (8 sample positions)

Material of Construction: MS powder-coated

Digestion Unit (cms) : 33 x 28 x 75 (8 sample positions) (L x W x H) (Inclusive of clamp)

Plot No 39 1st Floor Lavanya Aracade Jayabheri Enclave Gachibowli Hyderabad-500032.Telangana
Mob : +91 99082 46685 Tel : 040 - 35172288/89