



Digital Overhead Stirrer

FEATURES:

1. Independent LED display for RPM & Timer 99:59 (HH:MM)
2. User Friendly
3. Brushless DC motor
4. Liquid cannot reach internal components due to the front safeguard
5. Audible alarm for error conditions
6. Effective stirring provides clockwise & anti-clockwise stirring modes
7. Prominently used in chemical synthesis, pharmaceuticals, physical & chemical analysis, petrochemical, cosmetics, health products, foods, biotechnology & so on



TECHNICAL SPECIFICATIONS:

Model No.	6100.RSolv	6050.RSolv
T-On Time	5-999 Sec	5-999 Sec
T-off Time	0, 5-999 Sec	0, 5-999 Sec
Speed Range	50-3000 RPM	50-1500 RPM
Wifi Control	NO	NO
Temperature Measurement	NO	NO
Max. Viscosity	6000 mPass	10000 mPass
Max. Torque	65 N.CM	80 N.CM
Time Range	1 min-99.59 min	1 min-99.59 min
Max. Stir Capacity (H2O)	10-40 L	40-60 L
Stirring Shaft Extension Feature	Yes	Yes
Stirring Modes	Yes	Yes
Motor Type	BLDC	BLDC
Operating Voltage	24 VDC	24 VDC
Supply Voltage	100-240 VAC, 50/60Hz	100-240 VAC, 50/60Hz
Power	130	130
Dimension	90X240X220 mm	90X240X220 mm
Net Weight	12 Kgs.	12 Kgs.



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.