



RIKI GLOBAL
performance through innovation

Vortex Shaker

(CYCLO MIXER)



DESCRIPTION:

Test tube shaker, vibration type, variable speed, vortex genie short rapid eccentric stroke and motor driven. Neoprene cap imparts vigorous agitation to tube or small flasks. Shaking starts when vessel is pressed against cup, with ON/OFF switch and indicator, vacuum rubber fitted at the bottom to prevent movement.



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.