



RIKI GLOBAL
performance through innovation



Magnetic stirrer with Hot Plate

RS-PSI introduced world's best Magnetic-stirrers. This group of magnetic stirrers is similar to the conventional magnetic stirrers but have additional hot plate to achieve heating of solutions to be stirred. The heating is controlled by energy regulator, which ensures regulation at any desired percentage of the capacity of heater.

Technical Specifications:

Maximum Speed RPM: 1800

Stirring Capacity Ltrs: 2 liter

Heating Capacity Watts 500

Dimension (W x D x H) mm 150 x150 x

170 Stirring Paddle ØxL mm: 25 x 6