



# Orbital Shaking Incubator

## DESCRIPTIONS:

Our company also takes pride in providing a wide range of orbital shaker These are designed & fabricated to suit various is standard of orbital shaking motion applicable in growing field of medical agricultural industrial Research Laboratories & Others.

## FEATURES:

1. Double-walled inner chamber made of S.S.304 & outer M. S. Powder coated
2. PUF insulation between two walled
3. Heavy angle frame structure from all sides
4. Shaking assembly electric pulley mechanism
5. Platform adjustable stainless steel platform
6. Flask Holder: Spring-type clamp arrangement with easy interchangeability facility Heavy duty Motor (DC Motor)
7. RPM 50 to 250 RPM controlled by Regulator
8. Stroke 25 to 30 mm Stroke displacement
9. Temp. Controlled by digital Temp. Controller cum indicator with accuracy of  $\pm 0.5$  °C
10. Temp. Range 10 °C to 60 °C in the cooling Model & in the heating Model 5 °C above ambient to 60 °C
11. Cooling by CFC-free hermetically sealed compressor
12. Heating by 'U' Type tubular S. S. Heater Illumination light to view
13. RPM Indicator in digital Display
14. Digital timer with Audio visual Alarm

## SIZE AVAILABLE FOR ORBITAL INCUBATOR:

No.of Flask X Size in ML	Rating Watts
9 X250 ml	500
16 X250 ml	1000
25 X250 ml	1500
9 X250 ml	500
16 X250 ml	1000
25 X250 ml	1500



Plot No 39, 1st Floor, Lavanya Arcade,  
Jayabheri Enclave, Gachibowli,  
Hyderabad-500032.Telangana.