



Viscometer

DESCRIPTION:

RSolv Viscometer is a scientific instrument used to measure the viscosity of a fluid. Viscosity is a measure of a fluid's resistance to deformation or flow, and it describes how easily a fluid can flow. In simple terms, it is a measure of a fluid's thickness or stickiness.

FEATURES:

1. Full color, touch screen display
2. Accurate and reliable operation
3. Direct readout of all measurement parameters
4. Viscosity, Temperature, % Torque, Spindle/Speed
5. Auto Range Showing
6. Time Function for measurement
7. Sound alarm at under 20% Torque
8. Linear calibration



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.