


Specifications of RSolv - Capella HPLC Phases

COLUMN DIMENSIONS:

Columns are available in lengths of 50, 75, 100, 150 and 250 mm and with i.d.s of 2.0, 4.0 and 4.6 mm. Semi-preparative and Preparative HPLC columns can be supplied with lengths of 50, 100, 150 and 250 mm and with i.d.s of 7.8, 10.0, 21.2, 30.0 and 50.0 mm



Phase	Particle	Pore Size	Surface Area	Carbon Load	Endcapped	USP	Phase Code
C18	3, 5, 0	100	325	19	Yes	L1	01
C8	3, 5, 0	100	325	11	Yes	L7	02
C6	3, 5, 0	100	325	8	Yes	L15	03
C4	3, 5, 0	100	325	6	No	L26	04
Phenyl	3, 5, 0	100	325	12	Yes	L11	05
PFP	3, 5, 0	100	325	9	Yes	L43	06
CN	3, 5, 0	100	325	6	No	L10	07
AMINO	3, 5, 0	100	325	4	No	L8	08
DIOL	3, 5, 0	100	325	4	No	L20	09
Silica	3, 5, 0	100	325	n/a	n/a	L3	10



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.