



Bagmixer

DESCRIPTIONS:

Avoiding all risks of cross-contamination, RS-BagMixer 400W is an easy-to-use a powerful lab blender. Adapted to all kind of applications and with a gurantee of optimal bacterial extraction, it is a great tool for lab analyses.

1. QUICK HOMOGENIZING OF SAMPLES

- Optimal bacterial extraction
- Only 60 seconds to homogenize a sample

2. EASY CLEANING

- Fully opening door (270°)
- Security drip tray
- All stainless steel

3. USER-FRIENDLY

- Compact
- Adjustable blending time (30 -210 secs. or open-running)
- Q-Tight® closing: no more leaks upwards

SPECIFICATIONS:

Brand	Interscience
Model	RS-BagMixer 400 W
Capacity	50-400 mL
Material	Stainless Steel
Blending speed	Fixed speed 8 strokes/s
Dimensions	26 x 39 x 29 cm
Weight	16.5 kg
Power	110-220 V 50-60Hz
Blending time	30,60,90,120,150,180,210 s
Compact and ergonomic	Yes
Leak free Q tight system	Yes
Window door	Yes
Security drip tray	Yes
Door	270 deg opening door
Door weight	950 g



Bagmixer 400W



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.