




Digital Colony Counter

DESCRIPTIONS:

Digital colony counter is designed for quick and accurate counting of bacterial and mould colonies in petri dishes. Our digital colony counter is a feature packed ft easy to use, tins is an indispensable bench top tool for the busy microbiologist. It is designed for rapid and accurate counting of bacterial and mould colonies.

FEATURES:

1. 100 mm Dia Magnifying Glass
2. Audible Confirmation of Each Count
3. Uniform Glair Free illumination
4. Digital Readout from 0-9999
5. Auto Marker Pen
6. Wolffhuegel Glass Grid with Focusing Facility



SPECIFICATIONS:

Display	4 Digit, 9999 Max. Count
Dish Size	110 mm
Magnification	x 1.7
Dimension	l 274 X b 320 X h 167 mm (approx.)
Weight	3 Kgs. (approx.)
Packing	23 cm X 28 cm X 20 cm Corrugated Carton
Power	230 V - 10 % AC, 50 Hz, 40 W
Accessories	Marking Pen - 1 No. Magnifier Lens - 1 No.



Plot No 39, 1st Floor, Lavanya Arcade,
Jayabheri Enclave, Gachibowli,
Hyderabad-500032.Telangana.