



RIKI GLOBAL
performance through innovation

KJELDAHL DIGESTER



FEATURES

- Soft touch tactile membrane keyboard with LED indication

SPECIFICATIONS

- Heating block solid: Aluminum
- No. of samples: 8 nos
- Glass Flask: 250 ml
- Temperature : 10° C above ambient to 400° C over temp
- Sample Size: 0.25 - 5 gms
- Protection Available Temp: Resolution 0.1° C sensor PT 100
- Timer: User definable up to 99 Hrs 59 Min
- Display: 7-segment LED keyboard soft touch membrane
- Power: 230V Ac, 50 Hz wattage 1500 Watts (8 sample positions)
- Material of Construction: MS powder-coated
- Digestion Unit (cms) : 33 x 28 x 75 (8 sample positions) (L x W x H)
(Inclusive of clamp)