



RS-SC-2 AUTOMATIC SEED COUNTING MACHINE

DESCRIPTION

Automatic seed Counter can automatic count many kinds of grains by setting number. In agriculture it is mainly used in grain count, exp. Rice, wheat, sorghum, corn, and vegetables seed and so on.



PRODUCT FEATURES

- Two modes: Free count and bulk count
- The counting speed is adjustable
- High accuracy: Sensitivity levels for different sizes of particles
- Automatic deceleration: when the counting number is close to the setting number, the counting speed slows down automatically
- Automatic motor Protection: the motor will stop working once it detects there are no particles on the plate for 60 seconds
- User-friendly interface: LCD screen, mode selection, counting time, level, and number of particles displayed clearly
- Wide application: suitable for particles includes but is not limited to crop seeds, vegetable seeds, feed, diamond, and tablets

TECHNICAL PARAMETER

Count Range	0.7-12mm
Count Accuracy	±2‰
Count Speed	≥1000 /3min (Sesame)
Counting Capacity	1~99999
Bulk Counting	1-99999
Noise	<70db
Power	AC110V /60Hz Or AC 220V /50Hz
Continuous working time	≥ 5Hr
Size of Machine	285*210*176mm
Operation Environment	-10°C- 50°C, <85% RH